



多沃生物

Dowobio Biotechnology Co., Ltd



## 山羊抗小鼠 HRP 二抗

产品货号	产品名称	包装
DW1068	山羊抗小鼠 HRP 二抗	100ul

### Background:

The secondary antibodies are available conjugated to enzyme, biotin or fluorophore for use in a variety of antibody-based applications including Western Blot, ImmunoHistoChemistry, ImmunoFluorescence, Flow Cytometry and ELISA. We offer high quality secondary antibodies from goat, rabbit and donkey sources for your each application. Serum adsorbed secondary antibodies are also available and are recommended for use with immunoglobulin-rich samples.

### Reactivity :

Mouse

### Applications :

WB,IHC-p,IF(paraffin section),ELISA

### Formulation :

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

### Source :

Goat

### Dilution :

WB 1:5000-20000, ELISA 1:10000-20000, IHC 1:500-1000

### Purification :

The antibody was affinity-purified from goat antiserum by affinity- chromatography using epitope-specific immunogen.

### Concentration :

4 mg/ml

### Storage Stability :

-20°C/1 year

### 本产品文章中的写法:

英文: Dowobio (Shanghai, China)

中文: 上海多沃生物科技有限公司 Dowobio, 上海