




Mono-2-ethylhexyl phthalate-[d4]

产品货号	产品名称	包装
DW3089	Mono-2-ethylhexyl phthalate-[d4]	2.5mg

产品介绍:

Description	Mono-2-ethylhexyl phthalate-[d4] is the labelled analogue of Monoethylhexyl phthalic acid (MEHP). MEHP is a metabolite of Bis(2-ethylhexyl)phthalate (DEHP).
Synonyms	Mono-2-ethylhexyl phthalate-D4; Mono(2-ethylhexyl) Phthalate-d4; 2-Ethylhexyl Hydrogen Phthalate-d4; Mono(2-ethylhexyl) Ester; 1,2-(Benzene-d4)dicarboxylic Acid Mono(2-ethylhexyl) Ester; MEHP-d4; Mono(ethylhexyl) Phthalate-d4
IUPAC 名称	2,3,4,5-tetradeuterio-6-(2-ethylhexoxycarbonyl)benzoic acid
CAS	1276197-22-8
Canonical SMILES	<chem>CCCCC(CC)COC(=O)C1=CC=CC=C1C(=O)O</chem>
Molecular Weight:	282.37
Molecular Formula:	 $C_{16}H_{18}D_4O_4$
Purity	95% by HPLC; 98% atom D
Solubility	Soluble in Chloroform (Slightly), Methanol (Slightly)
Appearance	Colourless to Light Beige Oil to Semi-Solid
Storage	Store at 2-8°C

本产品文章中的写法:

英文: Dowobio (Shanghai, China)

中文: 上海多沃生物科技有限公司 Dowobio, 上海